



Gauze, white petroleum impregnated, sterile, overwrap strip, 3" x 36"



- Also known as Vaseline gauze, petroleum gauze
- Similar to Xeroform gauze

- Applied directly to wound with the intent of the gauze becoming adherent
- Trim as needed
- Can be used in conjunction with Silvadene, applying Silvadene first and the petroleum gauze over top as a primary dressing

INDICATIONS FOR USE

- Can be used for partial thickness injuries – moist, blanching wounds

Acticoat 40 cm x 40 cm sheet

- Uses microsilver technology to maintain broad-spectrum antimicrobial coverage
- Thin, low-adherent, multi-layered dressing
- Outer layers are coated with silver, while inner layer is absorbent



- Do not use on patients with a known sensitivity to silver
- Remove for MRI imaging
- May cause transient discoloring of surrounding skin
- Can be cut to size
- Three-five day wear time
- Needs to remain moist (NOT soaked) with sterile water (NOT sodium chloride) to remain active
- Either side can be applied to wound

INDICATIONS FOR USE

- Can be used for partial thickness injuries – moist, blanching wounds
- Can be used for full thickness injuries prior to the operating room
- Ideal for patients who have difficulty performing dressing changes



Silver sulfadiazine 1% cream, 400 mg jar

- Also written as SSD or Silvadene
- Readily available topical antimicrobial that provides coverage for many gram-negative and gram-positive bacteria (except for pseudomonas)

- Should be changed at least BID, but daily at minimum
- Use caution in 3rd trimester of pregnancy and first two months of life related to hyperbilirubinemia
- Can cause neutropenia, discontinue if this develops
- Avoid use in patients with a sulfa allergies
- The cream is expected to turn gray when oxidized. This does not affect its efficacy.

INDICATIONS FOR USE

- Any burn wounds
- Any open wounds



Burn dressing, sterile, 18" x 18" or 36" x 36"

- Comes in 2 sizes
- Multilayered
- Essentially a really large gauze



- Can be applied over Silvadene cream or Acticoat and secured with bandage roll or tubular dressing retainer

18" x 18"

- Adult upper extremity
- Adult Calf
- Pediatric upper and/or lower extremity
- Pediatric chest/back
- Head/Scalp

36" x 36"

- Adult chest/back
- Adult lower extremity

Bandage roll, sterile, 4.5" x 4.1 yds, 6 ply

- Also known as Kerlix or gauze roll



- ALWAYS wrap from most distal point of extremity to just beyond the wound
- ALWAYS wrap digits separately
- Caution on wrapping dressings too tightly

Elastic dressing retainer, tubular, roll, SIZES 1, 3, 5, 8, and 10

- Also known as stretchy, Stretch Net, Spandage, Elastinet, Flexinet



SIZE 1

Best for fingers

SIZE 3

Pediatric arm

Infant Leg

SIZE 5

Adult forearm

Pediatric leg

SIZE 8

Adult arm

Head

Adult calf

Pediatric trunk

SIZE 10

Adult leg

Adult trunk

- Secondary dressing to secure primary dressing in place

Solution bowl, plastic, 7 qt.

- Large clean bowl to be used for wound care



Solution bowl, plastic, sterile, 16 oz

- Small sterile bowl to be used for wound care



Towel drape, sterile, 18" and 26"



Suture removal kit



Glove, pair, sterile, surgical, non-latex, powder-free, SIZES 7, 7.5, and 8

